

## ALLURI (778XLSBN / TÜV Ready)





14-15	1350	59-62
<small>Speed (km/h)</small>	<small>Power (Watt)</small>	<small>Range (km)</small>

Rating: Not Rated Yet

[Ask a question about this product](#)

#### Description

The Shoprider 778XLSBN is a large scooter with a weight capacity of 225 kg. Beautifully crafted and sturdily built, this is yet another addition to SHOPRIDER's unique line of scooters. The Delta Style Tiller is sleekly molded with a textured vinyl finish. The floor of the scooter is designed with room enough for the larger rider as is the 20" deluxe seat which is standard.

#### \*For the German market:

**This so called TÜV-ready model, needs to be presented to a TÜV institute of your choice, before getting the type approval (BE = Betriebserlaubnis) which is needed to drive the vehicle in road traffic. For further or more detailed information, please contact your Shoprider Dealer near you.**