

VETAS (888IX)





9-10	860	26-30
<small>Speed (km/h)</small>	<small>Power (Watt)</small>	<small>Range (km)</small>

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

This scooter offers a modern streamline design with many updated features including an easy to use tiller dash, full LED lighting with super bright headlight, gas adjustable tiller and 265 mm flat free low profile tyres. The 888iX is also fitted with comfort ride suspension and a high fully rotating, lockable seat. The reliability of these impressive scooters by SHOPRIDER™ is further enhanced by a powerful motor, rear wheel drive and electromagnetic brake.

The 888iX has the ability to travel wide open footpaths. Or take you right through your local shopping centre.