

RUNNER 2 (888SLB / TÜV Ready)





10-11	1280	36-40
Speed (km/h)	Power (Watt)	Range (km)

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

The SHOPRIDER™ 888SLB is compact in design and lightweight in style, and guaranteed to provide the proven safety and reliable features of a Shoprider scooter.

With its functional design, the 888SLB is very comfortable and will provide you with a high level of personal independence and freedom.

The reliability of the 888SLB is enhanced by the rear wheel drive and electromagnetic brake.

Features front, rear and hazard lights, a padded high back seat, finger or thumb controls and a basket on the front tiller. The front and rear suspensions insure a comfortable ride.

***For the German market:**

This so called TÜV-ready model, needs to be presented to a TÜV institute of your choice, before getting the type approval (BE = Betriebserlaubnis) which is needed to drive the vehicle in road traffic.

For further or more detailed information, please contact your Shoprider Dealer near you.