

## AGILIA (888WA)



## Wheel Chair : AGILIA (888WA)

---



7-8	580	14-16
<small>Speed (km/h)</small>	<small>Power (Watt)</small>	<small>Range (km)</small>

Rating: Not Rated Yet

[Ask a question about this product](#)

### Description

The Shoprider 888WA Power Wheelchair combines style and functionality. It comes with puncture proof 6"/10" tires, 5A Off-Board charger and a sloped shroud.

It also offers an adjustable reclining seat, adjustable footrest, and 3 way adjustable armrests.

Perfect for hallways - only 55.88cm wide, it has a small turning radius of 650 mm which makes it one of the most manoeuvrable and popular power chairs.