

## AVIDI (888WNLLHD)



5-6	600x2	26-28
<small>Speed (km/h)</small>	<small>Power (Watt)</small>	<small>Range (km)</small>

Rating: Not Rated Yet  
[Ask a question about this product](#)

#### Description

The Shoprider 888WNLLHD has been ergonomically engineered to provide the robust combination of style, comfort and extra weight capacity (205kg).

It has the flexibility of a multi adjustable 510mm captains seat i.e. seat cushion height, seat cushion angle, backrest angle & forward or aft adjustment on seat base, when combined with the powerful smooth motors, mid wheel manoeuvrability this electric powerchair has everything it needs.

The Shoprider 888WNLLHD's smooth operation is further enhanced by the easy steer & operate joystick controller, 360mm puncture proof tyres and front and rear elastometers.

The Shoprider 888WNLLHD is suited to operation on a variety of terrains and hard surfaces, both inside and outside, the 50Amp batteries fitted also ensure it can travel long distances between charges.