

MAXI CAB (889XLSBN-CABIN)





14-14,5	1350	40-42
<small>Speed (km/h)</small>	<small>Power (Watt)</small>	<small>Range (km)</small>

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

Introducing SHOPRIDER's top of line scooter. The ultra exclusive 889XLSBN-CABIN will satisfy those with discriminating taste for only the very best.

Expecting inclement weather? No problem! The cabin scooter's stylish canopy enclosure provides comfortable, all season transportation.

Sturdy, key locking and removable doors provide security as well as convertibility.

The 889XLSBN-CABIN also comes with a fully operational lighting system, full suspension for a comfortable ride, even a windscreen wiper/washer system with reservoir.

Comes with a high quality heating system as an option.