

ALIZE (779XLS)



Large : ALIZE (779XLS)

---



|              |           |            |
|--------------|-----------|------------|
| 18           | 1350      | 50-54      |
| Speed (km/h) | Power (W) | Range (km) |

Rating: Not Rated Yet  
[Ask a question about this product](#)

Description

**brand-new chopper style design**

The 779XLS model is designed in a chopper style, with bare pipes and visible shock absorbers. At the same time, it provides the comfort of a large mobility scooter with an upgraded front suspension system. This nice looking 3-wheel scooter is very manoeuvrable and easy to drive. It has a weight capacity of 185 kg combined with a high back contoured sliding/rotating/height adjustable seat with headrest.