

## BRIZA (889ASN)





11-12	750	30-32
Weight (kg)	Power (W)	Range (km)

Rating: Not Rated Yet

[Ask a question about this product](#)

#### Description

### Stylish 4-WHEELER

This cool looking new model comes loaded with all the features you expect from a Shoprider scooter including Gas Adjustable Tiller, Powerful 1.5 hp motor, Soft Ride Suspension and Airflow Comfort Seat. The Briza is further enhanced by large wheels for improved ground clearance, Steering Tiller with easy-to-use controls, rear wheel drive and electromagnetic brake.

### Features

- Gas adjustable steering tiller
- comfortable handle with finger operated controls
- Front + Rear (Adjustable) suspension
- Full light package: Twin Headlights (LED), Rear lights, Hazard lights
- Rear lights– USB power outlet