

## CITY STAR & CITY STAR 2 (GK10)





5-6	780	33-35
<small>Speed (km/h)</small>	<small>Power (W)</small>	<small>Range (km)</small>

Rating: Not Rated Yet

[Ask a question about this product](#)

#### Description

The GK10 is the largest of the SHOPRIDER™ portable mobility scooters. This is a scooter that is suited to a person that wants all the autonomy and comfort of a full sized scooter but requires the convenience of an easily detachable one

This scooter is fitted with 36 amp batteries that can be charged with the scooter or the battery pack can be removed and charged separately while the scooter remains in storage.

The Gk 10 detachable scooter separates into four modules that will easily fit into the boot of a car so the machine is transportable. The heaviest module is only 21 kilograms

To separate the modules there are quick release levers so no tools are required.