

SWIFTY (GK4)



Compact : SWIFTY (GK4)



5-7	360	10-12
<small>Speed (km/h)</small>	<small>Power (Watt)</small>	<small>Range (km)</small>

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

After the SL74, the GK4 is the most compact of the SHOPRIDER™ portable mobility scooters. Perfect for travel and easy for storage.

This scooter is fitted with 12-14 amp batteries that can be charged with the scooter or the battery pack can be removed and charged separately while the scooter remains in storage.

The GK4 model includes solid tyres, a comfortable swivel seat with arm rests, and a basket.