

## RUNNER 2 (888SLB / TÜV Ready)





10-11	1280	36-40
<small>Speed (km/h)</small>	<small>Power (Watt)</small>	<small>Range (km)</small>

Rating: Not Rated Yet

[Ask a question about this product](#)

#### Description

The SHOPRIDER™ 888SLB is compact in design and lightweight in style, and guaranteed to provide the proven safety and reliable features of a Shoprider scooter.

With its functional design, the 888SLB is very comfortable and will provide you with a high level of personal independence and freedom.

The reliability of the 888SLB is enhanced by the rear wheel drive and electromagnetic brake.

Features front, rear and hazard lights, a padded high back seat, finger or thumb controls and a basket on the front tiller. The front and rear suspensions insure a comfortable ride.

**\*For the German market:**

**This so called TÜV-ready model, needs to be presented to a TÜV institute of your choice, before getting the type approval (BE = Betriebserlaubnis) which is needed to drive the vehicle in road traffic.**

**For further or more detailed information, please contact your Shoprider Dealer near you.**